

Antibodies Product Datasheet

Product Name:	CHCHD4 Antibody
Catalog Number:	C43062
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CHCHD4
Concentration:	1 mg/ml
Gene Name:	CHCHD4 MIA40
Protein Name:	CHCHD4
Other Name:	Mitochondrial intermembrane space import and assembly protein 40 (Coiled-coil-helix-coiled-coil-helix domain-containing protein 4)
Human Swiss Prot No.:	Q8N4Q1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.03% sodium azide.
Specificity:	This antibody detects endogenous levels of CHCHD4 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year