

Antibodies Product Datasheet

Product Name:	CHD6 Antibody
Catalog Number:	C43066
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	CHD6 CHD5 KIAA1335 RIGB
Protein Name:	Chromodomain-helicase-DNA-binding protein 6 (CHD-6) (EC 3.6.4.12) (ATP-dependent helicase CHD6) (Radiation-induced gene B protein)
Other Name:	
Human Swiss Prot No.:	Q8TD26
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CHD6 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year