

Antibodies Product Datasheet

Product Name:	CHD9 Antibody
Catalog Number:	C43067
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	CHD9 KIAA0308 KISH2 PRIC320 AD-013 x0008
Protein Name:	Chromodomain-helicase-DNA-binding protein 9 (CHD-9) (EC 3.6.4.12) (ATP-dependent helicase CHD9) (Chromatin-related mesenchymal modulator) (CReMM) (Chromatin-remodeling factor CHROM1) (Kismet homolog 2) (PPAR-alpha-interacting complex protein 320 kDa) (Per
Other Name:	
Human Swiss Prot No.:	Q3L8U1
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CHD9 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year