

Antibodies Product Datasheet

Product Name:	CHIP Antibody
Catalog Number:	C43071
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 160-240
Concentration:	1 mg/ml
Gene Name:	STUB1 CHIP PP1131
Protein Name:	E3 ubiquitin-protein ligase CHIP (EC 6.3.2.-) (Antigen NY-CO-7) (CLL-associated antigen KW-8) (Carboxy terminus of Hsp70-interacting protein) (STIP1 homology and U box-containing protein 1)
Other Name:	
Human Swiss Prot No.:	Q9UNE7
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CHIP Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year