

Antibodies Product Datasheet

Product Name:	CHM4B Antibody
Catalog Number:	C43075
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CHM4B
Concentration:	1 mg/ml
Gene Name:	CHMP4B C20orf178 SHAX1
Protein Name:	Charged multivesicular body protein 4b (Chromatin-modifying protein 4b) (CHMP4b) (SNF7 homolog associated with Alix 1) (SNF7-2) (hSnf7-2) (Vacuolar protein sorting-associated protein 32-2) (Vps32-2) (hVps32-2)
Other Name:	
Human Swiss Prot No.:	Q9H444
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CHM4B
Storage and Stability:	-20°C/1 year