

Antibodies Product Datasheet

Product Name:	CHST5 Antibody
Catalog Number:	C43078
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	CHST5
Protein Name:	Carbohydrate sulfotransferase 5 (EC 2.8.2.-) (Galactose/N-acetylglucosamine/N-acetylglucosamine 6-O-sulfotransferase 4-alpha) (GST4-alpha) (Intestinal N-acetylglucosamine-6-O-sulfotransferase) (I-GlcNAc6ST) (Intestinal GlcNAc-6-sulfotransferase) (hIGN6ST)
Other Name:	
Human Swiss Prot No.:	Q9GZS9
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CHST5 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year