

## Antibodies Product Datasheet

|                               |                                                                                                                                                |
|-------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | CHSTB Antibody                                                                                                                                 |
| <b>Catalog Number:</b>        | C43080                                                                                                                                         |
| <b>Reactivity:</b>            | Human, Mouse, Rat                                                                                                                              |
| <b>Applications:</b>          | WB, ELISA                                                                                                                                      |
| <b>Dilutions:</b>             | WB 1:500-2000 ELISA 1:5000-20000                                                                                                               |
| <b>Modification:</b>          | Unmodified/Total                                                                                                                               |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                                              |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                          |
| <b>Immunogen:</b>             | Synthesized peptide derived from part region of human protein                                                                                  |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                                        |
| <b>Gene Name:</b>             | CHST11                                                                                                                                         |
| <b>Protein Name:</b>          | Carbohydrate sulfotransferase 11 (EC 2.8.2.5) (Chondroitin 4-O-sulfotransferase 1) (Chondroitin 4-sulfotransferase 1) (C4S-1) (C4ST-1) (C4ST1) |
| <b>Other Name:</b>            |                                                                                                                                                |
| <b>Human Swiss Prot No.:</b>  | Q9NPF2                                                                                                                                         |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.                                                                                 |
| <b>Specificity:</b>           | CHSTB Polyclonal Antibody detects endogenous levels of protein.                                                                                |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                                   |