

Antibodies Product Datasheet

Product Name:	CHSTC Antibody
Catalog Number:	C43081
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	CHST12 UNQ500/PRO1017
Protein Name:	Carbohydrate sulfotransferase 12 (EC 2.8.2.5) (Chondroitin 4-O-sulfotransferase 2) (Chondroitin 4-sulfotransferase 2) (C4ST-2) (C4ST2) (Sulfotransferase Hlo)
Other Name:	
Human Swiss Prot No.:	Q9NRB3
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CHSTC Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year