

Antibodies Product Datasheet

Product Name:	CHSTF Antibody
Catalog Number:	C43083
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 440-520
Concentration:	1 mg/ml
Gene Name:	CHST15 BRAG GALNAC4S6ST KIAA0598
Protein Name:	Carbohydrate sulfotransferase 15 (EC 2.8.2.33) (B-cell RAG-associated gene protein) (hBRAG) (N-acetylgalactosamine 4-sulfate 6-O-sulfotransferase) (GalNAc4S-6ST)
Other Name:	
Human Swiss Prot No.:	Q7LFX5
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CHSTF Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year