

Antibodies Product Datasheet

Product Name:	CL054 Antibody
Catalog Number:	C43122
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CL054
Concentration:	1 mg/ml
Gene Name:	C12orf54 HSD-29 HSD-30 HSD29 HSD30
Protein Name:	Uncharacterized protein C12orf54
Other Name:	
Human Swiss Prot No.:	Q6X4T0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CL054
Storage and Stability:	-20°C/1 year