

Antibodies Product Datasheet

Product Name:	CL12A Antibody
Catalog Number:	C43124
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	CLEC12A CLL1 DCAL2 MICL
Protein Name:	C-type lectin domain family 12 member A (C-type lectin-like molecule 1) (CLL-1) (Dendritic cell-associated lectin 2) (DCAL-2) (Myeloid inhibitory C-type lectin-like receptor) (MICL)
Other Name:	
Human Swiss Prot No.:	Q5QGZ9
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CL12A Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year