

Antibodies Product Datasheet

Product Name:	CLASR Antibody
Catalog Number:	C43130
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	CLASRP SFRS16 SWAP2 UNQ2428/PRO4988
Protein Name:	CLK4-associating serine/arginine rich protein (Splicing factor, arginine/serine-rich 16) (Suppressor of white-apricot homolog 2)
Other Name:	
Human Swiss Prot No.:	Q8N2M8
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CLASR Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year