

Antibodies Product Datasheet

Product Name:	CLC2D Antibody
Catalog Number:	C43134
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 40-120
Concentration:	1 mg/ml
Gene Name:	CLEC2D CLAX LLT1 OCIL
Protein Name:	C-type lectin domain family 2 member D (Lectin-like NK cell receptor) (Lectin-like transcript 1) (LLT-1) (Osteoclast inhibitory lectin)
Other Name:	
Human Swiss Prot No.:	Q9UHP7
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CLC2D Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year