

Antibodies Product Datasheet

Product Name:	CLC4D Antibody
Catalog Number:	C43135
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	CLEC4D CLECSF8 MCL
Protein Name:	C-type lectin domain family 4 member D (C-type lectin superfamily member 8) (C-type lectin-like receptor 6) (CLEC-6)
Other Name:	
Human Swiss Prot No.:	Q8WXI8
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CLC4D Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year