

Antibodies Product Datasheet

Product Name:	CLC7A Antibody
Catalog Number:	C43138
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	CLEC7A BGR CLECSF12 DECTIN1 UNQ539/PRO1082
Protein Name:	C-type lectin domain family 7 member A (Beta-glucan receptor) (C-type lectin superfamily member 12) (Dendritic cell-associated C-type lectin 1) (DC-associated C-type lectin 1) (Dectin-1)
Other Name:	
Human Swiss Prot No.:	Q9BXN2
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CLC7A Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year