

Antibodies Product Datasheet

Product Name:	CLC9A Antibody
Catalog Number:	C43139
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CLC9A
Concentration:	1 mg/ml
Gene Name:	CLEC9A UNQ9341/PRO34046
Protein Name:	C-type lectin domain family 9 member A
Other Name:	
Human Swiss Prot No.:	Q6UXN8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CLC9A
Storage and Stability:	-20°C/1 year