

Antibodies Product Datasheet

Product Name:	CLCN5 Antibody
Catalog Number:	C43144
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	CLCN5 CLCK2
Protein Name:	H(+)/Cl(-) exchange transporter 5 (Chloride channel protein 5) (ClC-5) (Chloride transporter ClC-5)
Other Name:	
Human Swiss Prot No.:	P51795
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CLCN5 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year