

Antibodies Product Datasheet

Product Name:	CLVS1 Antibody
Catalog Number:	C43161
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CLVS1
Concentration:	1 mg/ml
Gene Name:	CLVS1 CRALBPL RLBP1L1
Protein Name:	CLVS1
Other Name:	Clavesin-1 (Cellular retinaldehyde-binding protein-like) (Retinaldehyde-binding protein 1-like 1)
Human Swiss Prot No.:	Q8IUQ0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.49% sodium azide.
Specificity:	This antibody detects endogenous levels of CLVS1 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year