

Antibodies Product Datasheet

Product Name:	CM035 Antibody
Catalog Number:	C43163
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CM035
Concentration:	1 mg/ml
Gene Name:	C13orf35
Protein Name:	Putative uncharacterized protein C13orf35
Other Name:	
Human Swiss Prot No.:	Q6ZP68
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CM035
Storage and Stability:	-20°C/1 year