

Antibodies Product Datasheet

Product Name:	CMIP Antibody
Catalog Number:	C43165
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CMIP
Concentration:	1 mg/ml
Gene Name:	CMIP KIAA1694 TCMIP
Protein Name:	CMIP
Other Name:	C-Maf-inducing protein (c-Mip) (Truncated c-Maf-inducing protein) (Tc-Mip)
Human Swiss Prot No.:	Q8IY22
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.345% sodium azide.
Specificity:	This antibody detects endogenous levels of CMIP at Human/Mouse
Storage and Stability:	-20°C/1 year