

Antibodies Product Datasheet

Product Name:	CNNM4 Antibody
Catalog Number:	C43189
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CNNM4
Concentration:	1 mg/ml
Gene Name:	CNNM4 ACDP4 KIAA1592
Protein Name:	CNNM4
Other Name:	Metal transporter CNNM4 (Ancient conserved domain-containing protein 4) (Cyclin-M4)
Human Swiss Prot No.:	Q6P4Q7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.358% sodium azide.
Specificity:	This antibody detects endogenous levels of CNNM4 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year