

Antibodies Product Datasheet

Product Name:	CNOTA Antibody
Catalog Number:	C43190
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CNOTA
Concentration:	1 mg/ml
Gene Name:	CNOT10
Protein Name:	CCR4-NOT transcription complex subunit 10
Other Name:	
Human Swiss Prot No.:	Q9H9A5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CNOTA
Storage and Stability:	-20°C/1 year