

Antibodies Product Datasheet

Product Name:	COASY Antibody
Catalog Number:	C43209
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 260-340
Concentration:	1 mg/ml
Gene Name:	COASY PSEC0106
Protein Name:	Bifunctional coenzyme A synthase (CoA synthase) (NBP) (POV-2) [Includes: Phosphopantetheine adenylyltransferase (EC 2.7.7.3) (Dephospho-CoA pyrophosphorylase) (Pantetheine-phosphate adenylyltransferase) (PPAT); Dephospho-CoA kinase (DPCK) (EC 2.7.1.24) (D
Other Name:	
Human Swiss Prot No.:	Q13057
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	COASY Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year