

Antibodies Product Datasheet

Product Name:	COE3 Antibody
Catalog Number:	C43211
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 50-130
Concentration:	1 mg/ml
Gene Name:	EBF3 COE3
Protein Name:	Transcription factor COE3 (Early B-cell factor 3) (EBF-3) (Olf-1/EBF-like 2) (O/E-2) (OE-2)
Other Name:	
Human Swiss Prot No.:	Q9H4W6
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	COE3 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year