

Antibodies Product Datasheet

Product Name:	COG7 Antibody
Catalog Number:	C43213
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human COG7
Concentration:	1 mg/ml
Gene Name:	COG7 UNQ3082/PRO10013
Protein Name:	COG7
Other Name:	Conserved oligomeric Golgi complex subunit 7 (COG complex subunit 7) (Component of oligomeric Golgi complex 7)
Human Swiss Prot No.:	P83436
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.332% sodium azide.
Specificity:	This antibody detects endogenous levels of COG7 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year