

Antibodies Product Datasheet

Product Name:	COL12 Antibody
Catalog Number:	C43214
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	COLEC12 CLP1 NSR2 SCARA4 SRCL
Protein Name:	Collectin-12 (Collectin placenta protein 1) (CL-P1) (hCL-P1) (Nurse cell scavenger receptor 2) (Scavenger receptor class A member 4) (Scavenger receptor with C-type lectin)
Other Name:	
Human Swiss Prot No.:	Q5KU26
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	COL12 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year