

Antibodies Product Datasheet

Product Name:	COLA1 Antibody
Catalog Number:	C43215
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human COLA1
Concentration:	1 mg/ml
Gene Name:	COL21A1 COL1AL FP633
Protein Name:	COLA1
Other Name:	Collagen alpha-1(XI) chain
Human Swiss Prot No.:	Q96P44
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.244% sodium azide.
Specificity:	This antibody detects endogenous levels of COLA1 at Human
Storage and Stability:	-20°C/1 year