

Antibodies Product Datasheet

Product Name:	COLGALT1 Antibody
Catalog Number:	C43216
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human COLGALT1
Concentration:	1 mg/ml
Gene Name:	GLT25D1 PSEC0241
Protein Name:	COLGALT1
Other Name:	Procollagen galactosyltransferase 1 (EC 2.4.1.50) (Glycosyltransferase 25 family member 1) (Hydroxylysine galactosyltransferase 1)
Human Swiss Prot No.:	Q8NBJ5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.04% sodium azide.
Specificity:	This antibody detects endogenous levels of COLGALT1 at Human/Mouse
Storage and Stability:	-20°C/1 year