

Antibodies Product Datasheet

Product Name:	COX7R Antibody
Catalog Number:	C43229
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 10-90
Concentration:	1 mg/ml
Gene Name:	COX7A2L COX7AR COX7RP
Protein Name:	Cytochrome c oxidase subunit 7A-related protein, mitochondrial (COX7a-related protein) (Cytochrome c oxidase subunit VIIa-related protein) (EB1)
Other Name:	
Human Swiss Prot No.:	O14548
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	COX7R Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year