

Antibodies Product Datasheet

Product Name:	COX7S Antibody
Catalog Number:	C43230
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human COX7S
Concentration:	1 mg/ml
Gene Name:	COX7A2P2 COX7A3 COX7AL2 COX7AP2
Protein Name:	Putative cytochrome c oxidase subunit 7A3, mitochondrial (Cytochrome c oxidase subunit VIIa 3)
Other Name:	
Human Swiss Prot No.:	O60397
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human COX7S
Storage and Stability:	-20°C/1 year