

Antibodies Product Datasheet

Product Name:	CP062 Antibody
Catalog Number:	C43236
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CP062
Concentration:	1 mg/ml
Gene Name:	C16orf62 101F10.2
Protein Name:	CP062
Other Name:	UPF0505 protein C16orf62 (Esophageal cancer-associated protein)
Human Swiss Prot No.:	Q7Z3J2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.261% sodium azide.
Specificity:	This antibody detects endogenous levels of CP062 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year