

Antibodies Product Datasheet

Product Name:	CP27B Antibody
Catalog Number:	C43240
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	CYP27B1 CYP1ALPHA CYP27B
Protein Name:	25-hydroxyvitamin D-1 alpha hydroxylase, mitochondrial (EC 1.14.13.13) (25-OHD-1 alpha-hydroxylase) (25-hydroxyvitamin D(3) 1-alpha-hydroxylase) (VD3 1A hydroxylase) (Calcidiol 1-monooxygenase) (Cytochrome P450 subfamily XXVIIIB polypeptide 1) (Cytochrome
Other Name:	
Human Swiss Prot No.:	O15528
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CP27B Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year