

Antibodies Product Datasheet

Product Name:	CP2C9 Antibody
Catalog Number:	C43241
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	CYP2C9 CYP2C10
Protein Name:	Cytochrome P450 2C9 (EC 1.14.13.-) ((R)-limonene 6-monooxygenase) (EC 1.14.13.80) ((S)-limonene 6-monooxygenase) (EC 1.14.13.48) ((S)-limonene 7-monooxygenase) (EC 1.14.13.49) (CYP11C9) (Cytochrome P-450MP) (Cytochrome P450 MP-4) (Cytochrome P450 MP-8) (C
Other Name:	
Human Swiss Prot No.:	P11712
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CP2C9 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year