

Antibodies Product Datasheet

Product Name:	CP4F3 Antibody
Catalog Number:	C43247
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 400-480
Concentration:	1 mg/ml
Gene Name:	CYP4F3 LTB4H
Protein Name:	Leukotriene-B(4) omega-hydroxylase 2 (EC 1.14.13.30) (CYPIVF3) (Cytochrome P450 4F3) (Cytochrome P450-LTB-omega) (Leukotriene-B(4) 20-monooxygenase 2)
Other Name:	
Human Swiss Prot No.:	Q08477
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CP4F3 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year