

Antibodies Product Datasheet

Product Name:	CPT2 Antibody
Catalog Number:	C43260
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CPT2
Concentration:	1 mg/ml
Gene Name:	CPT2 CPT1
Protein Name:	Carnitine O-palmitoyltransferase 2, mitochondrial (EC 2.3.1.21) (Carnitine palmitoyltransferase II) (CPT II)
Other Name:	
Human Swiss Prot No.:	P23786
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CPT2
Storage and Stability:	-20°C/1 year