

Antibodies Product Datasheet

Product Name:	CQ10A Antibody
Catalog Number:	C43270
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CQ10A
Concentration:	1 mg/ml
Gene Name:	COQ10A UNQ6192/PRO20219
Protein Name:	Coenzyme Q-binding protein COQ10 homolog A, mitochondrial
Other Name:	
Human Swiss Prot No.:	Q96MF6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CQ10A
Storage and Stability:	-20°C/1 year