

Antibodies Product Datasheet

Product Name:	CR3L4 Antibody
Catalog Number:	C43277
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 200-280
Concentration:	1 mg/ml
Gene Name:	CREB3L4 AIBZIP CREB4 JAL
Protein Name:	Cyclic AMP-responsive element-binding protein 3-like protein 4 (cAMP-responsive element-binding protein 3-like protein 4) (Androgen-induced basic leucine zipper protein) (AIBZIP) (Attaching to CRE-like 1) (ATCE1) (Cyclic AMP-responsive element-binding pro
Other Name:	
Human Swiss Prot No.:	Q8TEY5
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CR3L4 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year