

Antibodies Product Datasheet

Product Name:	CRDL2 Antibody
Catalog Number:	C43279
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 330-410
Concentration:	1 mg/ml
Gene Name:	CHRD2 BNF1 CHL2 UNQ765/PRO1557
Protein Name:	Chordin-like protein 2 (Breast tumor novel factor 1) (BNF-1) (Chordin-related protein 2)
Other Name:	
Human Swiss Prot No.:	Q6WN34
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CRDL2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year