

Antibodies Product Datasheet

Product Name:	CREL2 Antibody
Catalog Number:	C43283
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CREL2
Concentration:	1 mg/ml
Gene Name:	CRELD2 UNQ185/PRO211
Protein Name:	CREL2
Other Name:	Cysteine-rich with EGF-like domain protein 2
Human Swiss Prot No.:	Q6UXH1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.32% sodium azide.
Specificity:	This antibody detects endogenous levels of CREL2 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year