

Antibodies Product Datasheet

Product Name:	CRHBP Antibody
Catalog Number:	C43287
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 70-150
Concentration:	1 mg/ml
Gene Name:	CRHBP CRFBP
Protein Name:	Corticotropin-releasing factor-binding protein (CRF-BP) (CRF-binding protein) (Corticotropin-releasing hormone-binding protein) (CRH-BP)
Other Name:	
Human Swiss Prot No.:	P24387
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CRHBP Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year