

Antibodies Product Datasheet

Product Name:	CRIS2 Antibody
Catalog Number:	C43288
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 30-110
Concentration:	1 mg/ml
Gene Name:	CRISP2 GAPDL5 TPX1
Protein Name:	Cysteine-rich secretory protein 2 (CRISP-2) (Cancer/testis antigen 36) (CT36) (Testis-specific protein TPX-1)
Other Name:	
Human Swiss Prot No.:	P16562
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CRIS2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year