

Antibodies Product Datasheet

Product Name:	CRLF3 Antibody
Catalog Number:	C43292
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 270-350
Concentration:	1 mg/ml
Gene Name:	CRLF3 CREME9 CYTOR4 P48
Protein Name:	Cytokine receptor-like factor 3 (Cytokine receptor-like molecule 9) (CREME-9) (Cytokine receptor-related protein 4) (Type I cytokine receptor-like factor p48)
Other Name:	
Human Swiss Prot No.:	Q8IUI8
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CRLF3 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year