

Antibodies Product Datasheet

Product Name:	CRNN Antibody
Catalog Number:	C43294
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 380-460
Concentration:	1 mg/ml
Gene Name:	CRNN C1orf10 DRC1 PDRC1 SEP53
Protein Name:	Cornulin (53 kDa putative calcium-binding protein) (53 kDa squamous epithelial-induced stress protein) (58 kDa heat shock protein) (Squamous epithelial heat shock protein 53) (Tumor-related protein)
Other Name:	
Human Swiss Prot No.:	Q9UBG3
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CRNN Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year