

Antibodies Product Datasheet

Product Name:	CRYAA Antibody
Catalog Number:	C43299
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 1-80
Concentration:	1 mg/ml
Gene Name:	CRYAA CRYA1 HSPB4
Protein Name:	Alpha-crystallin A chain (Heat shock protein beta-4) (HspB4) [Cleaved into: Alpha-crystallin A chain, short form]
Other Name:	
Human Swiss Prot No.:	P02489
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CRYAA Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year