

Antibodies Product Datasheet

Product Name:	CSN7A Antibody
Catalog Number:	C43310
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	COPS7A CSN7A DERP10
Protein Name:	COP9 signalosome complex subunit 7a (SGN7a) (Signalosome subunit 7a) (Dermal papilla-derived protein 10) (JAB1-containing signalosome subunit 7a)
Other Name:	
Human Swiss Prot No.:	Q9UBW8
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CSN7A Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year