

Antibodies Product Datasheet

Product Name:	CSPG4 Antibody
Catalog Number:	C43311
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 1910-1990
Concentration:	1 mg/ml
Gene Name:	CSPG4 MCSP
Protein Name:	Chondroitin sulfate proteoglycan 4 (Chondroitin sulfate proteoglycan NG2) (Melanoma chondroitin sulfate proteoglycan) (Melanoma-associated chondroitin sulfate proteoglycan)
Other Name:	
Human Swiss Prot No.:	Q6UVK1
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CSPG4 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year