

Antibodies Product Datasheet

Product Name:	CSTN2 Antibody
Catalog Number:	C43314
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CSTN2
Concentration:	1 mg/ml
Gene Name:	CLSTN2 CS2
Protein Name:	CSTN2
Other Name:	Calsyntenin-2 (Alcadein-gamma) (Alc-gamma)
Human Swiss Prot No.:	Q9H4D0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.235% sodium azide.
Specificity:	This antibody detects endogenous levels of CSTN2 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year