

Antibodies Product Datasheet

Product Name:	CSTN3 Antibody
Catalog Number:	C43315
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CSTN3
Concentration:	1 mg/ml
Gene Name:	CLSTN3 CS3 KIAA0726
Protein Name:	CSTN3
Other Name:	Calsyntenin-3 (Alcadein-beta) (Alc-beta)
Human Swiss Prot No.:	Q9BQT9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.237% sodium azide.
Specificity:	This antibody detects endogenous levels of CSTN3 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year