

Antibodies Product Datasheet

Product Name:	CTBL1 Antibody
Catalog Number:	C43324
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CTBL1
Concentration:	1 mg/ml
Gene Name:	CTNNBL1 C20orf33 PP8304
Protein Name:	CTBL1
Other Name:	Beta-catenin-like protein 1 (Nuclear-associated protein) (NAP) (Testis development protein NYD-SP19)
Human Swiss Prot No.:	Q8WYA6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.107% sodium azide.
Specificity:	This antibody detects endogenous levels of CTBL1 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year