

Antibodies Product Datasheet

Product Name:	CTDP1 Antibody
Catalog Number:	C43326
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CTDP1
Concentration:	1 mg/ml
Gene Name:	CTDP1 FCP1
Protein Name:	CTDP1
Other Name:	RNA polymerase II subunit A C-terminal domain phosphatase (EC 3.1.3.16) (TFIIF-associating CTD phosphatase)
Human Swiss Prot No.:	Q9Y5B0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.255% sodium azide.
Specificity:	This antibody detects endogenous levels of CTDP1 at Human/Mouse
Storage and Stability:	-20°C/1 year